Platelet Poor Plasma

REF 5343022 Platelet Poor Plasma 5 x 2 mL

Symbols key / Symbolschlüssel / символы
(if applicable / wenn anwendbar / если используются)

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PRODUCT DESCRIPTION**INTENDEND USE**

Platelet Poor Plasma can be used as a negative control for the Lupus Anticoagulant Test and the TECHNOCLOT® LA Screen and Confirm dRVVT assay.

COMPOSITION

The Platelet Poor Plasma is prepared from normal plasma. The platelet content of the plasma is reduced to a level where it does not influence the Lupus tests by careful handling of the plasma and repeated centrifugation.

MATERIAL REQUIRED (not supplied with the kit)

- Pipette
- Distilled water
- Required Lupus Tests

REF	5343005	Lupus Anticoagulant Test	6 tests
REF	5279025	CaCl ₂ solution 50 mmol/L	100 mL

REF	5343012	TECHNOCLOT® LA Screen	5x2 mL
REF	5343016	TECHNOCLOT® LA Confirm	5x1 mL

- Required positive control

REF	5343010	Lupus Inhibitor Plasma	2x1 mL
REF	5343019	Lupus Inhibitor Plasma	5x1 mL
REF	5343021	Lupus Inhibitor Plasma Low	5x1 mL

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Lupus Inhibitor Plasma are tested and found negative for Hb_sAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT* (Ceveron**)	2...8°C	-20°C
1 day	3 days	1 month

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

* = room temperature

** = in the Ceveron® alpha in the respective control area in the sample tray

TEST PROCEDURE**CEVERON**

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

PREPARATION AND PERFORMANCE

- Open the vial carefully and reconstitute the contents in 2 mL of distilled water by carefully rotating the vial (avoid frothing).
- Allow the reconstituted Plasma to stand for 10 minutes at room temperature before use. For standardisation a reconstitution time of 30 min is recommended.
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

ANALYSES RESULTS

See instructions for Lupus Anticoagulant Test and the TECHNOCLOT® LA Screen and Confirm dRVVT assay. Please refer to the values for LA screen, LA confirm and LA Ratio stated in the lot specific batch table.

STANDARDISATION

No international calibrator is available for the standardization of Lupus Inhibitor tests. The Lupus Inhibitor documentation is based on studies testing plasma samples from normal subjects and plasma samples with lupus inhibitor. A house reference batch is established in order to avoid batch-to-batch variation.

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