Trinity Biotech

TriniLIA™ Auto Dimer

REF T3102	1 Kit
------------------	--------------

INTENDED USE

TriniLIA Auto Dimer is an immunoturbidimetric assay used for the quantitative determination of the fibrin degradation products that contain D-Dimer in human plasma on automated instruments.

SUMMARY AND PRINCIPLE

D-dimer containing moieties are formed by plasmin degradation of factor XIIIa cross-linked fibrin. Elevated levels of D-Dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC).¹⁻³ Laboratory measurements of fibrin degradation products, including D-Dimer, have significance in the use of these conditions.

TriniLIA Auto Dimer is a turbidimetric assay that utilizes antibody coated latex particles. In the presence of D-Dimer, the particles aggregate and turbidity increases. The increase in scattered light is proportional to the amount of D-Dimer in the sample. The latex particles are coated with a monoclonal antibody that reacts with fibrin D-Dimer or fragment D of fibrin. The antibody has no cross reactivity with fibrinogen.⁴ This allows for the determination of D-Dimer in human plasma.

REAGENT

For *in vitro* diagnostic use.

REAGENT DESCRIPTION

The reagents are lot-specific. Lots are not interchangeable.

TriniLIA Auto Dimer Latex Reagent, 1 x 13.5 ml, T3102A

Latex particles coated with anti-D-Dimer monoclonal antibody MA-803, suspended in HEPES buffer pH 8.5, containing stabilizers, detergent and sodium azide.

TriniLIA Auto Dimer Low Standard, 1 x 22 ml, T3102B

HEPES buffer, pH 7.0, containing stabilizers, detergent and sodium azide.

TriniLIA Auto Dimer High Standard, 2 x 1 ml, T3102C

Lympholized human plasma immunodepleted of D-Dimer.

TriniLIA Auto Dimer High Standard, 2 x 1 ml, T3102D

Lympholized human plasma enriched with D-Dimer.

TriniLIA Auto Dimer Saline: 2 x 8 ml, T3000A

Anticoagulant, pH 7.3, containing sodium azide.

REAGENT PREPARATION

Latex Reagent, Reaction Buffer and Auto Dimer Saline are ready to use.

The latex may sediment during storage. Mix thoroughly before use.

Auto Dimer Low Standard

TriniLIA Auto Dimer Saline. The reagent is not used.

Other Instruments: Reconstitute with 1.0 ml of Auto Dimer Saline and agitate gently for 5 minutes to completely dissolve contents.

Auto Dimer High Standard

TriniLIA Auto Dimer Saline. Reconstitute with 1.4 ml of Auto Dimer saline and agitate gently for 5 minutes to completely dissolve content.

Note: This procedure will give the D-Dimer concentration of approximately 6.40 µgFEU/ml (2.40 mg/L). See lot specific document for the assay value.

Other Instruments: Reconstitute with 1.0 ml of Auto Dimer Saline and agitate for 5 minutes to completely dissolve content.

Note: This procedure will give the D-Dimer concentration of approximately 6.40 µgFEU/ml (320µg/L D-Dimer). See lot specific document for the assay value.

ADDITIONAL MATERIALS REQUIRED

Automatic analyzer

D-Pipettes (Low and High concentration),

Pipettes for 125 µl and 1000 µl

Plastic test tubes

MATERIALS AVAILABLE

TriniCHECK D-Dimer Low

TriniCHECK D-Dimer High

INSTRUMENTS

Trinity Biotech Instruments:

Applications/Method Adaptations for individual analyzers are available upon request, please contact your local representative.

The number of tests obtained with specific instruments may vary.

STORAGE AND STABILITY

Unopened reagents are stable until the expiration date provided on the vial label when stored at 2-8°C.

Latex Reagent and Reaction Buffer

Store at 2-8°C and use within 4 weeks from opening.

Instrument on-board storage:

Thrombolyzer application: Stable for 1 week when stored on-board.

Hitachi 911/902 application: Stable for 4 weeks when stored on-board.

Auto Dimer Saline

Store at 2-8°C and use within 4 weeks from opening.

After reconstitution:

Auto Dimer Low Standard and Auto Dimer High Standard Stable for 10 hours at 20-25°C.

SPECIMEN COLLECTION AND PREPARATION

Nine volumes of blood are to be collected in one volume of 3.2% (0.109M) sodium citrate. Immediately after blood collection, samples are centrifuged at 1500 x g for 15 minutes. Please refer to the manual for detailed methods. However no test can be performed unless the specimen is centrifuged at 1500 x g for 15 minutes. For further instructions regarding specimen collection and storage.⁵

PROCEDURE

WARNINGS AND PRECAUTIONS

TriniLIA Auto Dimer is of human origin. Each donor unit of source plasma used in these products has been tested and found negative for Hepatitis B antigens, HIV I and II antibodies, Hepatitis C antibodies, syphilis antibodies and H.T.L.V. II antibody. TriniLIA Auto Dimer is not a blood product. It does not contain any animal products. It does not contain any complete assurance that products derived from human blood will not transmit infectious disease. As with all materials of human origin this product should be handled as a potentially infectious agent. All wastes containing biological material should be properly labelled and stored separately from other wastes. Dispose of all waste materials according to prescribed international, national and local regulations.

The Latex Reagent, Reaction Buffer and Auto Dimer Saline contain sodium azide, which may react with lead or copper plumbing to form lead or copper azides. Materials discarded into a sink should be flushed with a large volume of water to prevent azide build-up.

The test should be used in conjunction with clinical observations and results of other laboratory tests.

TEST PROCEDURE

Preparation of Standard Curve

Prepare the standard dilutions as described below. The standard dilutions are stable for 10 hours at 20-25°C.

Users must construct a standard curve each time a new kit lot is used, or every 6 months, and if control Plasma is assayed out of range.

Each standard curve shows are shown in the instrument-specific application sheets.

Thrombolyzer application: The dilutions are made by the instrument using Auto Dimer Saline solution as diluent. The following dilutions should be programmed:

Programmed Dilution	Dilution Factor	Low standard (µg/l)	Mix with
A	1	11.3	-
B	2	125	125 µl Std Dilution A
C	4	125	125 µl Std Dilution B
D	8	125	125 µl Std Dilution C
E	16	125	125 µl Std Dilution D
G	32	125	125 µl Std Dilution E
G	0	125	-

Use the value for Auto Dimer High Standard given in the lot-specific document to determine the exact D-Dimer concentration in each standard dilution.

The D-Dimer concentration is calculated by dividing the lot specific assigned value by the "dilution factor" corresponding to each dilution.

For example: The lot specific assigned value of the high standard is 6.66 µgFEU/ml (332µg/L) then:

the concentration corresponding to dilution B is: 6.66 µgFEU/ml (332µg/L)÷2=3.33 µgFEU/ml (166µg/L).

- the concentration corresponding to dilution C is: 6.66 µgFEU/ml (332µg/L)÷4=1.66 µgFEU/ml (833 µg/L).

QUALITY CONTROL

It is recommended that the control plasmas TriniCHECK D-Dimer Low and TriniCHECK D-Dimer High are assayed at regular intervals in order to ensure consistent assay results. If the control plasma results deviates from the D-Dimer concentration given in the lot-specific document, a new standard curve should be constructed.

RESULTS

INTERPRETATION OF RESULTS

Conversion of results from mg/L to µg/L: 1 mg/L = 1000 µg/L.

The results may be reported in purified D-Dimer or Fibrinogen Equivalent Units.

All information in this document is given in µg FEU/ml and µg/L purified DD in brackets.

If necessary to convert ng/ml or µg/L D-Dimer unit to FEU use the following formula:

1000 ng/ml purified D-Dimer = 2 µgFEU/ml

The following conversion chart can be utilized to convert D-Dimer results. This table is to be used only as a general comparison.

ng/ml or µg/L D-Dimer	ngFEU/ml or µgFEU/L	µg/ml or mg/L D-Dimer	µgFEU/ml
100	200	0.100	0.20
300	600	0.3	0.6
1000	2000	1	2
2000	4000	2	4

EXPECTED RESULTS

In a study of 80 (Hitachi 911) and 30 (Hitachi 902) normal individuals, 95% of the values were below 0.22 µgFEU/ml (110 µg/L).

In a study of 30 normal individuals on Thrombolyzer, 95% of the values were below 0.30 µgFEU/ml (150 µg/L).

Elevated levels are found in patients with confirmed Deep Venous Thrombosis (DVT), Pulmonary Embolism (PE), and Trauma. D-Dimer levels rise during pregnancy and high levels are associated with complications.⁵

The concentration of D-Dimer in any given specimen may differ from the concentration determined using D-Dimer assays from different manufacturers. Consequently, each laboratory must determine reference intervals for their individual test populations, reagents, and instruments.

LIMITATIONS

Presence of rheumatoid arthritis factor may result in false-positive results influencing no interference.

On MA-803, suspended in HEPES buffer (6.7 µg/L), or triglycerides (17 µg/L). Highly lipemic samples should be diluted in Auto Dimer Saline and re-assayed.

On TriniCHECK, no interference is found from bilirubin (<0.27 µg/L), hemoglobin (<6.7 g/L), or triglycerides (3.3 g/L). Highly lipemic samples should be diluted in Auto Dimer Saline and re-assayed.

Human anti-mouse antibody testing was not performed with this product. Results from patients with heterophilic antibody should be interpreted with caution since this kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values.

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TriniLIA Auto Dimer, when used on different analyzers, correlated as follows:

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