smart CRP Test

In Vitro Diagnostic test kit for determination of c-reactive protein in serum or plasma and whole blood with smart photometer.

Eurolyser Diagnostica GmbH
Bayernstraße 11a
5020 Salzburg/Austria
Tel. +43 (0) 662 43 21 00
Fax. +43 (0) 662 43 21 00 50
www.eurolyser.com

Order Information	Indication	Kit size
Order No.: ST0100	smart CRP-Test Kit	32 tests/pack
Order No.: ST1000	smart CRP-Control Kit	2 x 1 ml (low high range)
Order No.: ST1010	smart CRP-Calibrator Kit	1 x 1 ml

English:

Test kit preparation:
Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary:
C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

Methode-Measurement Range:
Immunturbidimetric Assay:
Serum/Plasma: 0,5-120mg/l for smart546 and 1,0-120mg/l for smart700/340 or smart700
Whole blood: 2,0-240mg/l for smart546 and 2,0-240mg/l for smart700/340 or smart700

Principle:
Kinetic determination of the concentration of CRP by photometric measurement at 546nm or 700nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

Test Kit:
ERS Cuvette pre-filled with 1000µl buffer reagent Glycine Buffer 170mM
ERS Cap pre-filled with 250µl latex reagent with rabbit anti human CRP antibodies (0,2%)

Storage Instructions and Stability of the Test Kit:
The reagents are stable until the indicated month of expiry , if stored at 2-8°C.

Warnings and Precautions:
The reagents contain sodium azide (0,95g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.. Take the necessary precautions for the use of laboratory reagents.

Waste Management:
please refer to local legal requirements.

Specimen Sample Material:
Freshly drawn serum is preferred.
Heparin plasma, EDTA plasma, edta whole blood and capillary blood from finger tip may also be used if limited test performance is accepted.

Reference Range:
adults < 5mg/l (<0,5 mg/dl) ; values > 3mg/l (> 0,3 mg/dl) can indicate a risk for coronary heart disease CHD.
It is recommended that each laboratory establish its own expected range.

Quality Control:
for internal quality control the smart CRP Control Calibrator Kit order no.:ST1000 should be used

Precision for Serum:
within run: N=20 Mean=100 mg/l cv=3,5 %
correlation with Cobas 6000
y=0,938x+0,32 r=0,993 (n=50) x=Cobas6000 serum y=smartCRP whole blood

References:
1. Osmond, A.P. et al, Proc. Natl. Acad. Sci. 74:739-743,1977
2. Pepys, M.B. Lancet.1:653-657,1981

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2.) Attention: Select whole blood or serum via the analyser menu

- For detailed information please refer to the user manual of the smart laboratory photometer
- 2.1) Sample type menu
 - 2.2) Hematokrit
 - 2.3) Starting measurement

3.) Sample preparation for whole blood (fingertip-or EDTA primary tube)

- 3.1a) Aspirate 5µl from fingertip OR ...
- 3.1b) Aspirate 5µl from EDTA capillary blood collection system OR ...
- 3.1c) Aspirate 5µl from EDTA primary tube
- 3.2) Dispense 5µl sample into liquid in ERS cuvette
- 3.3) Apply ERS cap
- 3.4) Place ERS cuvette into smart analyser.

ATTENTION!

Do not wipe the pipette tip after aspirating the sample!
Incubate tests 10 minutes at roomtemperature before use!

4.) Sample preparation for Serum or Plasma (primary tube serum , Li-Hep or EDTA Plasma)

- 4.1a) Aspirate 5µl serum from primary tube OR ...
- 4.1b) Aspirate 5µl plasma from centrifuged capillary blood collection system
- 4.2) Dispense 5µl sample into liquid in ERS cuvette
- 4.3) Apply ERS cap
- 4.4) Place ERS cuvette into smart analyser.

5.) Sample processing with smart photometer

Please follow the instructions in the users manual