**FOR INFORMATION USE ONLY**  
**Not to be used for performing the assay**  
**Refer to the insert accompanying kit**

## MDA® Verify® Reference Plasma

**REF** 252573 **20 x 1.2 ml**

### INTENDED USE

MDA® Verify® Reference Plasma is an assayed human plasma that has been lyophilized to maintain the integrity of the constituents. It is intended for use as a reference plasma in the quantitation of thrombosis and hemostasis proteins.<sup>1,2</sup>

### SUMMARY AND PRINCIPLE

Many procedures performed in the hemostasis laboratory require the use of a well-characterized reference plasma to evaluate levels of clotting factors, inhibitors, and fibrinolytic proteins in blood specimens. It is common practice to use freshly drawn pools of normal platelet-poor plasma from a relatively large number of healthy donors as this reference plasma.<sup>3</sup>

It is difficult to obtain and maintain a consistent and reliable pool. Using a large number of qualified donors (those who are healthy, non-smoking, under no medication, etc.) does not ensure a reference material with 100% activity for each hemostasis protein. Therefore, a reliable reference plasma containing quantitated levels of the major hemostatic proteins is a suitable replacement for a normal pool. MDA® Verify® Reference Plasma is a stabilized, lyophilized reference plasma useful for preparing reference curves for clotting and chromogenic assays.

The assigned activity of each hemostasis protein included in MDA® Verify® Reference Plasma was calculated based on the results of multiple assays. All factors were assayed against WHO and other reference preparations in order to assure the reliability and accuracy of assigned values.

### REAGENT

For *in vitro* diagnostic use.

### REAGENT DESCRIPTION

**Caution:** Handle as if capable of transmitting infectious agents. This product was prepared from human plasma which has been tested using FDA-licensed methods and found to be nonreactive for the presence of HIV-1 and HIV-2 antibodies, hepatitis B surface antigen (HBsAg), as well as hepatitis C virus (HCV). However, as no test method can offer complete assurance that infectious agents are absent, all specimens of human origin, including patient samples, should be handled as though they contain potentially infectious agents.

20 x 1.2 ml MDA® Verify® Reference Plasma (human) – Lyophilized, buffered plasma containing protein stabilizers and optimal levels of hemostasis proteins.

### REAGENT PREPARATION

Prior to removing the safety seal and the vial stopper, tap the vial on the counter to ensure that the lyophilized material is concentrated on the bottom of the vial. Reconstitute MDA® Verify® Reference Plasma with 1.2 ml of purified water, USP or equivalent. Replace the stopper and swirl the vial gently. Allow to reconstitute undisturbed at room temperature for 2-3 minutes. Mix by gently inverting the vial 5-6 times. Allow to stand at room temperature for a total of 20 minutes. Remix prior to use by gently inverting the vial 5-6 times. Store reconstituted plasma at 2-15°C. Discard after 8 hours when using with all assays except for Factor VIII. Discard after 4 hours when using with Factor VIII. While each laboratory should establish its own time limitations based on assay, instrumentation, environment, and usage conditions, these limits should not exceed those stated above which have been found to be valid under controlled conditions.

### ADDITIONAL MATERIALS REQUIRED

Hemostasis instrumentation Control plasmas  
 Purified water, USP or equivalent Assay reagents

### MATERIALS AVAILABLE

Multi-Channel Discrete Analyzer (MDA)  
 MDA Assay Reagents  
 MDA Control Plasmas

### INSTRUMENTS

Use MDA® Verify® Reference Plasma with hemostasis instrumentation.

Review the instrument operator manual for additional information regarding the following:

1. Installation and special requirements.
2. Operation principles, instructions, precautions, and hazards.
3. Manufacturer's specifications and performance capabilities.
4. Service and maintenance information.

### STORAGE AND STABILITY

Store MDA® Verify® Reference Plasma at 2-8°C when not in use. The expiration date printed on the label indicates the date beyond which the unopened product should not be used. After the foil package has been opened, the product is stable for 6 months if kept at 2-8°C.

### PROCEDURE

#### WARNINGS AND PRECAUTIONS

1. Do not pipet any of the materials by mouth. Do not smoke, eat, or drink in areas in which specimens or reagents are handled.
2. Use disposable gloves and handle all blood specimens used in the assay as though capable of transmitting infectious agents. Consult a physician immediately in the event that blood materials are ingested or come in contact with open lacerations, lesions, mucous membranes or other breaks in the skin.
3. Immediately clean up any spillage of specimens with a 1:10 dilution of 5% sodium hypochlorite. Dispose of the cleaning material by an acceptable method.
4. Treatment of blood products prior to disposal:
  - a. Autoclave for 60 minutes at 121°C.
  - b. Incinerate in disposable materials.
  - c. Mix liquid waste with 5% sodium hypochlorite solution so that the final concentration is approximately 0.5% sodium hypochlorite. Allow to stand at least 30 minutes before disposal.

#### TEST PROCEDURE

Refer to the instrument operator manual for assay instructions.

Gently mix MDA® Verify® Reference Plasma prior to placing it on the instrument to ensure homogeneity. Use MDA® Verify® Reference Plasma in the same manner as a pool of fresh normal plasma.

Use MDA® Verify® Reference Plasma to prepare reference curves for hemostasis proteins according to the assay instructions for each protein.

#### QUALITY CONTROL

1. The use of control plasmas (MDA® Verify® 1, 2, and 3) is recommended for monitoring calibration assays following established laboratory quality control procedures. CLSI recommends that controls be tested with each group of assays, with a change in personnel, and with every 40 samples.<sup>6</sup>
2. If control values are out of range, do not report patient results. Determine which part of the instrument/reagent/control system is not functioning properly and correct it. After corrective measures are implemented and documented following good laboratory practice, retest the controls. If the results of the retest are within range, patient samples can be tested and reported.

### RESULTS

#### INTERPRETATION OF RESULTS

The reference values assigned to this lot of MDA® Verify® Reference Plasma (listed in the enclosed value chart) are expressed as percent of normal activity. Because the activity of some of the hemostasis proteins is not 100%, a correction factor must be taken into account. Some instruments, such as the MDA, automatically calculate and report the corrected activity. As an example, the correction factor for a protein with the assigned value of 90% is 0.90. These instruments automatically multiply the assayed patient value by the correction factor.

#### Sample calculations

If a patient sample assays at 80% of normal for a protein that has an assigned activity of 90% with MDA® Verify® Reference Plasma, the value automatically reported by some instruments, such as the MDA, would be determined using the following formula:

$$\text{Assayed patient value} \times \text{Correction factor} = \text{Reported value} \\ 80\% \times 0.90 = 72\%$$

### LIMITATIONS

Variability within hemostasis assay results is known to occur and can be influenced by the type of instrumentation, reagents, and factor deficient plasmas used.<sup>4,5</sup>

MDA® Verify® Reference Plasma is not recommended for use in the preparation of heparin standard curves.

### PERFORMANCE CHARACTERISTICS

#### Methods used to obtain assay values

The assigned activity of each hemostasis protein included in MDA® Verify® Reference Plasma was calculated based on the results of multiple independent determinations. The values were generated using Trinity Biotech reagents and instrumentation.

All factors contained in MDA® Verify® Reference Plasma were assigned based on WHO and other reference preparations.

For technical assistance in the U.S.A., contact Trinity Biotech Customer Service at 1-800-325-3424. Outside the U.S.A., contact your local Trinity Biotech Representative.

MDA and Verify are registered trademarks of Trinity Biotech in the USA and other countries.

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