## MDA® Thromboquik™

**FOR INFORMATION USE ONLY**  
**Not to be used for performing the assay.**  
**Refer to the insert accompanying kit.**

REF 252561

### INTENDED USE

MDA® Thromboquik™ is a thrombin reagent (bovine) for use in the determination of the Thrombin Time (TT) in human plasma.

### SUMMARY AND PRINCIPLE

The Thrombin Time assay is based on the ability of thrombin to catalyze the polymerization of fibrinogen into a fibrin clot. The mechanism involves the cleavage of peptides (fibrinopeptides A and B) from the amino-terminal ends of the alpha and beta chains of fibrinogen by the proteolytic action of thrombin.<sup>1</sup> The resulting fibrin monomer molecules ultimately polymerize into a fibrin clot.

The Thrombin Time assay is a qualitative screening assay used to detect abnormalities in this phase of coagulation. The assay is performed by adding a known quantity of low concentration thrombin reagent to a plasma specimen *in vitro* and measuring the time required for clot formation to occur.

Abnormalities affecting this stage of coagulation include quantitative and qualitative alterations in fibrinogen, increased fibrinolytic activity causing variations in Fibrin(ogen) Degradation Products (FDP), and interferences with fibrinogen polymerization.<sup>2</sup> The Thrombin Time assay is also sensitive to heparin and other circulating antithrombins.<sup>3,5</sup>

### REAGENT

For *in vitro* diagnostic use.

#### REAGENT DESCRIPTION

300 tests  
10 x 2.5 ml MDA® Thromboquik™ – Lyophilized bovine thrombin; contains HEPES buffer, calcium chloride, and stabilizers.  
(30 tests each)

#### REAGENT PREPARATION

Reconstitute MDA® Thromboquik™ with 2.5 ml of purified water, USP or equivalent. Mix gently to ensure complete rehydration. Discard after 7 days. While each laboratory should establish its own time limitations based on instrumentation, environment, and usage conditions, these limits should not exceed those stated above which have been found to be valid under controlled conditions.

#### ADDITIONAL MATERIALS REQUIRED

Hemostasis instrumentation  
Purified water, USP or equivalent  
Control plasmas

#### MATERIALS AVAILABLE

Multi-Channel Discrete Analyzer (MDA)  
MDA® Verify® Control Plasmas  
TriniCHECK Fibrinogen Low Control

#### INSTRUMENTS

Use MDA® Thromboquik™ with hemostasis instrumentation. Review the instrument operator manual for additional information regarding the following:

1. Installation and special requirements.
2. Operation principles, instructions, precautions, and hazards.
3. Manufacturer's specifications and performance capabilities.
4. Service and maintenance information.

### STORAGE AND STABILITY

Store MDA® Thromboquik™ at 2-8°C when not in use. The expiration date printed on the label indicates the date beyond which the unopened product should not be used. Signs of deterioration are reflected in quality control results outside the established laboratory range. After the foil package has been opened, the product remains suitable for use 60 days if stored at 2-8°C. Store reconstituted MDA® Thromboquik™ at 2-15°C.

### SPECIMEN COLLECTION AND PREPARATION

**Collection:** No special preparation or fasting of the patient is necessary. Blood should be anticoagulated with sodium citrate dihydrate at a starting concentration of 3.2% (0.109 M). To assure accurate results, a ratio of nine parts blood to one part anticoagulant (9:1) should be used. Recommended procedures for the collection of diagnostic blood specimens for coagulation testing have been published by the Clinical Laboratory Standards Institute (CLSI).<sup>5</sup> Specimens that demonstrate visible hemolysis should not be used. The results obtained with plasma specimens that are icteric or lipemic should be interpreted with caution.

**Storage:** Centrifuge capped specimens at 2500 x g for 15 minutes. If the plasma is to be frozen, remove from cells, then freeze rapidly (-20°C or lower) and thaw rapidly (37°C) to prevent denaturation of fibrinogen. It is recommended that testing be completed within the time limits specified in the following table.<sup>6</sup>

Sample Storage Temperature	Test Within
22-24°C	2 hours
2-4°C	4 hours
-20°C	2 weeks
-70°C	6 months

### PROCEDURE

#### WARNINGS AND PRECAUTIONS

1. Do not pipet any of the materials by mouth. Do not smoke, eat, or drink in areas where specimens or reagents are handled.
2. Use disposable gloves and handle all specimens used in the assay cautiously as though capable of transmitting infectious agents. Consult a physician immediately in the event that specimens are ingested or come in contact with open lacerations, lesions, mucous membranes, or breaks in the skin.
3. Immediately clean up any spillage of specimens with a 1:10 dilution of 5% sodium hypochlorite. Dispose of the cleaning material by an acceptable method.
4. Dispose of all specimens and materials used to perform the test as if they contain infectious agents. Treatment prior to disposal:
  - a. Autoclave for 60 minutes at 121°C.
  - b. Incinerate disposable materials.
  - c. Mix liquid waste with 5% sodium hypochlorite solution so that the final concentration is approximately 0.5% sodium hypochlorite. Allow to stand 30 minutes before disposal.

#### TEST PROCEDURE

Refer to the instrument operator manual for assay instructions.

When MDA® Thromboquik™ is used on the MDA, the following processes are automatically performed. The patient sample (0.075 ml) is pipetted into a cuvette and warmed to 37 ± 1°C. The instrument then adds 0.075 ml of warmed MDA® Thromboquik™ to the sample and simultaneously begins timing for clot formation.

#### PROCEDURAL NOTES AND PRECAUTIONS

Some investigators recommend that heparinized plasma samples that have been frozen not be used in the Thrombin Time procedure.<sup>7</sup>

#### QUALITY CONTROL

1. The use of two control plasmas (MDA® Verify® 1 and MDA® Verify® 2, MDA® Verify® 3 or TriniCHECK Fibrinogen Low Control) is recommended for monitoring the Thrombin Time assay following established laboratory quality control procedures. CLSI recommends that controls be assayed at the initiation of testing, at least once each shift, or with each group of assays.<sup>6</sup>
2. If control values are out of range, do not report patient results. Determine which part of the instrument/reagent/control system is not functioning properly and correct it. After corrective measures are implemented and documented following good laboratory practice, retest the controls. If they are within range, patient samples can be tested and reported.

### RESULTS

The time recorded (in seconds) is the patient's thrombin time.

#### INTERPRETATION OF RESULTS

Although MDA® Thromboquik™ may be used with most manual or automated methods for clot detection, each method may detect a slightly different endpoint. Caution must be used when comparing results from different methods.

#### EXPECTED RESULTS

Clothing times are dependent upon numerous factors including temperature, water quality, pH, ionic strength, test system, specimen collection and preservation, and patient population. As a guide to the user, TT data were obtained in 90 normal adults to establish a reference interval. Based on these results, the reference interval range was determined to be 10-13 seconds. This is only a guide to the user; specific TT reference intervals should be established by each laboratory.<sup>5</sup>

### PERFORMANCE CHARACTERISTICS

Within-run precision studies on both normal and heparinized plasma and controls using MDA® Thromboquik™ on hemostasis instrumentation gave Coefficients of Variation (CVs) of less than 3%. Day-to-day precision studies on normal and abnormal controls using similar instrumentation gave CVs of less than 5%. MDA® Thromboquik™ was shown to be sensitive to fibrinogen levels as low as 75 mg/dl.

For technical assistance in the U.S.A., contact Trinity Biotech Customer Service at 1-800-325-3424. Outside the U.S.A., contact your local Trinity Biotech Representative.

MDA and Verify are registered trademarks of Trinity Biotech in the USA and other countries. Thromboquik is a trademark of Trinity Biotech.

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