Trinity Biotech

MDA® Auto-Dimer®

REF	T3103	122 Tests / 1 Kit
INTENDED USE		

MDA Auto-Dimer is an immunoturbidimetric assay used for the quantitative determination of the fibrin degradation products that contain D-dimer in human plasma.

SUMMARY AND PRINCIPLE

Quantitative, automated Latex Agglutination Test for Fibrin D-dimer.

D-dimer containing moieties are formed by plasmin degradation of factor XIIIa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC).¹⁻³ Laboratory measurements of fibrin degradation products, including D-dimer, have significance in the initial assessment of these conditions.

MDA Auto-Dimer is a turbidimetric assay that utilizes antibody coated latex particles. In the presence of D-dimer, the particles aggregate and turbidity increases. The increase in scattered light is proportional to the amount of D-dimer in the sample. The latex particles are coated with a monoclonal antibody that reacts with fibrin D-dimer or fragment D of fibrin. The antibody has no cross reactivity with fibrinogen.⁴ This allows for the determination of D-dimer in human plasma.

REAGENT

For *In Vitro* Diagnostic Use

REAGENT DESCRIPTION

The reagents are lot-specific. Kit components are not interchangeable between kit lots.

- D-dimer Latex Reagent** 1 x 13,5 ml latex particles coated with anti-D-dimer monoclonal antibody MA-8D3 suspended in HEPES buffer pH 8.5, containing stabilizers, detergent and sodium azide.
- D-dimer Reaction Buffer** 1 x 22 ml HEPES buffer, pH 7.0, containing stabilizers, detergent and sodium azide.
- D-dimer Saline Solution** 2 x 8 ml buffered saline, pH 7,3, containing sodium azide.
- D-dimer Low Standard** 4 x 1 ml lyophilized human plasma immunodepleted of D-dimer.
- D-dimer High Standard** 4 x 1 ml lyophilized human plasma enriched with D-dimer. Refer to Lot Specific Information sheet for assigned value.
- D-dimer Low Control** 4 x 1 ml lyophilized human plasma enriched with D-dimer. Refer to Lot Specific Information sheet for assigned value.

REAGENT PREPARATION

- D-dimer Latex Reagent** The latex may sediment during storage. Mix thoroughly before use. Transfer the volume of reagent required to a blue capped MDA vial for use on the MDA.
- D-dimer Reaction Buffer** Transfer the volume of reagent required to a yellow capped MDA vial for use on the MDA.
- D-dimer Saline solution** Transfer the volume of reagent required to a green capped MDA vial for use on the MDA. Use to dilute high patient samples and for reconstitution of Auto-Dimer Control and Standards.
- D-dimer Low standard** Reconstitute with 1.0 ml of Saline Solution and agitate gently for 5 minutes to completely dissolve contents. Transfer the volume of standard required to a green capped MDA vial for use on the MDA. Use to dilute the high calibrator for preparation of the calibration curve.
- D-dimer High standard** Reconstitute with 1.0 ml of Saline Solution and agitate gently for 5 minutes to completely dissolve contents. Transfer the volume of standard required to a Pink capped MDA vial for use on the MDA.
- D-dimer Low Control** Reconstitute with 1.0 ml of Saline Solution and agitate gently for 5 minutes to completely dissolve contents. Transfer the volume of control required to a Pink capped MDA vial for use on the MDA.

ADDITIONAL MATERIALS REQUIRED

Pipettes
Pipette tips
Test tubes

MATERIALS PROVIDED

D-dimer Latex Reagent
D-dimer Reaction Buffer
D-dimer Saline Solution
D-dimer Standards, Low and High
D-dimer Low Control
Bag of 26 empty MDA vials
Pre-printed barcoded labels for each component.
Low Control and High Standard Lot Specific Information sheet.

INSTRUMENTS

The reagents of the **MDA Auto-Dimer** kit are intended to be used on the MDA instrument.

For detailed instructions, refer to the instrument specific application, AP106.

STORAGE AND STABILITY

The unopened and unreconstituted reagents are stable until the expiration date stated on the box and vial labels when stored at 2-8°C .

- D-dimer Latex Reagent** Store at 2-8°C and use within 4 weeks from opening. Stable for 7 days when stored on-board.
- D-dimer Reaction Buffer** Store at 2-8°C and use within 4 weeks from opening. Stable for 7 days when stored on-board.
- D-dimer Saline solution** Store at 2-8°C and use within 4 weeks from opening.
- D-dimer Low standard** Stable for 10 hours at 20-25°C.
- D-dimer High standard** Stable for 10 hours at 20-25°C.
- D-dimer Low Control** Stable for 1 week at 2-8°C. Stable for 24 hours when stored on-board.

SPECIMEN COLLECTION AND PREPARATION

Nine volumes of blood are to be collected in one volume of 3.2% (0.109 M) sodium citrate. Immediately after blood collection, samples are centrifuged at 1500 x g for 15 minutes. Please refer to the most recent version of the CLSI document H21 for further instructions regarding specimen collection and storage.

PROCEDURE

WARNINGS AND PRECAUTIONS

The **MDA Auto-Dimer** Standards are of human origin. Each donor unit of source plasma used in these products has been tested and found negative for Hepatitis B antigens, HIV I and II antibodies, Hepatitis C antibodies, syphilis antibodies and H.T.L.V. I/II antibodies by FDA approved methods. However no test can offer complete assurance that products derived from human blood will not transmit infectious disease. As with all materials of human origin, this product should be handled as a potentially infectious agent. All wastes containing biological material should be properly labeled and stored separately from other wastes. Dispose of all waste materials according to prescribed international, national and local regulations. The Latex Reagent, Reaction Buffer and Saline Solution contain sodium azide, which may react with lead and copper plumbing to form highly explosive metal azides. Materials discarded into a sink should be flushed with a large volume of water to prevent azide build-up. The test should be used in conjunction with clinical observations and results of other laboratory tests.

PREPARATION OF STANDARD CURVE

Refer to the **MDA Auto-Dimer** Application, Reference AP106.

TEST PROCEDURE

Refer to the **MDA Auto-Dimer** Application, Reference AP106.

QUALITY CONTROL

It is recommended that the D-dimer Low Control provided in the **MDA Auto-Dimer** kit and the D-dimer High Control from **Trinity Biotech** are assayed at regular intervals in order to ensure consistent assay results. If the control plasma result deviates from the D-dimer concentration given in the Lot-Specific Information sheet a new standard curve should be constructed and controls re-run.

RESULTS

The results are reported in µg FEU/ml D-dimer.

EXPECTED RESULTS

Elevated levels are found in patients with confirmed deep venous thrombosis (DVT), pulmonary embolism, DIC, and trauma.¹⁻³ D-dimer levels rise during pregnancy and high levels are associated with complications.⁹ The concentration of D-dimer in any given specimen may differ from the concentration determined using D-dimer assays from different manufacturers. Consequently, each laboratory must determine reference intervals for their individual test populations, reagents, and instruments.

Please refer to the MDA Auto-Dimer application for expected values and performance characteristics for the MDA Instrument.

LIMITATIONS

Presence of rheumatoid arthritis factor may result in false-positive results (influence not quantified). No interference is found from bilirubin (<20 mg/dl) or hemoglobin (<1250 mg/dl g/L). Results from lipemic samples should be interpreted with caution. Highly lipemic samples should be diluted in Saline Solution and reassayed. Human anti-mouse antibody testing was not performed with this product. Results from patients with heterophilic antibody should be interpreted with caution since this kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values.

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REFERENCES / RIFERIMENTI / ΑΝΑΦΟΡΕΣ
<ol style="list-style-type: none">Declerck, P., et al. Fibrinolytic response and fibrin fragment D-dimer levels in patients with deep vein thrombosis. Thrombosis and Haemostasis 58, 1024-1029, 1987. Lindahl,T., et al. Clinical evaluation of a diagnostic strategy for deep venous thrombosis with exclusion by low plasma levels of fibrin degradation product D-dimer. Scand. J. Clin. Lab. Invest. 58, 307-316, 1998. Hansson, P.O., et al. Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? J.Intern. Med. 235, 143-151, 1994. Holvoet, P., et al. Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. Thrombosis and Haemostasis 61, 307-313, 1989. Ballegeer, V., et al. Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. Thrombosis and Haemostasis 58, 1030-1032, 1987. Dempfle, C-E., et al. The Fibrin Assay Comparison Trial (FACT). Evaluation of 23 Quantitative D-dimer assays as basis for the development of D-dimer calibrators. Thrombosis and Haemostasis 85, 671-8, 2001. Clinical and Laboratory Standards Institute (CLSI). Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline—Fifth Edition. CLSI document H21-A5 Vol. 28, No. 5, 2008.
ORDERING INFORMATION / INFORMAZIONI PER GLI ORDINI / ΠΛΗΡΟΦΟΡΙΕΣ ΠΑΡΑΓΕΛΙΑΣ

KIT	MDA® Auto-Dimer®	
Catalogue No.	T3103	
Kit Content	Item	Quantity
	D-dimer Latex reagent	1 x 13.5 ml
	D-dimer Reaction buffer	1 x 22 ml
	D-dimer Saline Solution	2 x 8 ml
	D-dimer Low Standard	4 x 1 ml
	D-dimer High Standard	4 x 1 ml
	D-dimer Low Control	4 x 1 ml
	MDA Vials	1 x 26
ADDITIONAL REAGENTS AVAILABLE		
Catalogue No.	Item	Quantity
1468	D-dimer High Control	5 x 1 ml
150712	Saline Solution	5 x 8 ml
252586	MDA Empty Reagent Vials	20 x 3 ml

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