CRYOcheck™

WEAK LUPUS POSITIVE CONTROL PLASMA DE CONTROLE LUPUS POSITIF FAIBLE

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Bibliography / Bibliographie

1. Conley CL, Hartmann RC. A hemorrhagic disorder caused by circulating anticoagulant in patients with disseminated lupus erythematosus. J Clin Invest 1952; 31:621.
2. Feinstein DI, Rapaport SI. Acquired inhibitors of blood coagulation. Progress in hemostasis and thrombosis 1972; 1:75-95.
3. Schleider MA, Nachman RI, Jaffe EA, Coleman M. A clinical study of the lupus anticoagulant. Blood 1976; 48:499-509.
4. Bowie WEJ, Thompson JH, Pascuzzi CA. Thrombosis in systemic lupus erythematosus despite circulating anticoagulants. J Lab Clin Med 1963; 62:416-430.
5. Nilsson IM, Astedt B, Hedner U, Berezin D. Intrauterine death and circulating anticoagulant, "Antithromboplastin". Acta Med Scand 1975; 197:153-159.
6. Thiagarajan P, Shapiro SS, DeMarco L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. J Clin Invest 1980; 66:397-405.
7. Triplett DA. Coagulation assays for the lupus anticoagulant: Review and critique of current methodology. Stroke 1992; (suppl. 1); 23:11.
8. Permpikul P, Rao LVM, Rapaport SI. Functional and binding studies of the roles of prothrombin and β_2 glycoprotein 1 in the expression of lupus anticoagulant activity. Blood 1994; 83:2878-2892.
9. Brandt JT, Triplett DA, Alving B, Scharrer I. Criteria for the diagnosis of lupus anticoagulants: an update. On behalf of the Subcommittee on Lupus Anticoagulants/ Antiphospholipid Antibody of the Scientific and Standardisation Committee of the ISTH. Thromb Haemost 1995; 74(4):1185-90.
10. Exner T, Rickard KA, Kronenberg H. A sensitive test demonstrating lupus anticoagulant and its behavioural patterns. Br J Haematol 1978; 40:143-151.
11. Thiagarajan P, Pengo V, Shapiro SS. The use of the dilute Russell's viper venom time for the diagnosis of lupus anticoagulants. Blood 1986; 68:869-884.
12. Alving BM, Baldwin PE, Richards RL, Jackson BJ. The dilute phospholipid APTT: A sensitive assay for verification of lupus anticoagulants. Thromb Haemost 1985; 54:709-712.
13. Triplett DA, Brandt JT, Kaczor D, Schaeffer J. Laboratory diagnosis of lupus inhibitors: A comparison of the tissue thromboplastin inhibition procedure with a new platelet neutralization procedure. Am J Clin Path 1983; 79:678-682.
14. Rauch J, Tannenbaum M, Tannenbaum H, Ramelson H, Cullis PR, Tilcock CP, et al. Human hybridoma lupus anticoagulants distinguish between lamellar and hexagonal phase lipid systems. J Biol Chem 1986; 261:9672-9677.
15. Biosafety in microbiological and biomedical laboratories 4th ed. Centers for Disease Control and Prevention / National Institutes of Health, 1999.
16. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
17. CLIA 1988 – Code of Federal Regulations, 42CFR493.1253, 1998.
18. Gilmer PR. Preanalytical variables in coagulation testing. In: Triplett, D. editor. Laboratory evaluation of coagulation. Illinois: College of American Pathologists; 1982. p. 1-8.

 European Authorized Representative (Regulatory affairs only)
Emergo Europe Molenstraat 15, 2513 BH The Hague, The Netherlands
Tel: (+31) 70 345 8570 Fax: (+31) 70 346 7299

Symbols used / Symboles utilisés

09.60.00016

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In vitro diagnostic medical device

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