

RVV-X

Description: Specific factor X activator from Russell's viper venom. RVV-X is a dimer of two peptide chains with a molecular weight of 60'000 g/mol each. Activation of factor X by RVV-X strictly depends on the presence of calcium ions.

Application: RVV-X is used in diagnostic procedures to quantitatively convert the zymogen factor X into factor Xa which can be determined by means of a clotting assay or photometrically, using a synthetic chromogenic substrate. RVV-X activator is used in testing of lupus anticoagulants.

Origin: *Vipera russelli* snake venom **MW:** approx. 120'000

Storage: May be used by the expiry date given on the label when stored unopened, protected from moisture, in the dark, 2°-8°C. Avoid contamination of the reagents by micro-organisms. Shipment of product does not require cooling during the time of transportation.

Assay :	<u>Chromogenic method</u>	<u>Clotting method</u>
	100 µl RVV-X (25 U/ml 25 mM CaCl ₂) 10 µl human citrated plasma => incubate for 75 s at 37 °C to activate factor X into factor Xa 790 µl 50 mM Tris-HCl, pH 8.4 100 µl Pefachrome®FXa (4 mM) => Determination of ΔOD/min at 405 nm	25 µl RVV-X (5 mU/ml 20 mM Tris- HCl, pH 7.2, 150 mM NaCl) 25 µl CaCl ₂ 25 mM 25 µl rabbit brain cephaline (0.2 mg/ ml) => incubate for 1 min at 37 °C 25 µl human citrated plasma

Unit definition: 1 unit (U) is the amount of RVV-X which generates one international unit of enzyme (factor Xa) from zymogen (factor X) [Stocker et al., 1986].

Stability after reconstitution:

+37°C	8 hours
+15 to +25°C	2 days
+2 to +8°C	1 week
-80°C	1 month


The reconstituted product can be frozen immediately and stored at -20°C.
Avoid refreezing.

References: Stocker K. Application of snake venom proteins in the diagnosis of hemostatic disorders. In: Medical Use of Snake Venom Proteins, Stocker K, ed. Boca Raton: CRC-Press 1990; 213-52.

Stocker K, Fischer H, Brogli M.
Determination of factor X activator in the venom of the saw-scaled viper (*Echis carinatus*),
Toxicon 1986; 24: 313-5.

Thiagarajan P, Pengo V, Shapiro SS. The use of dilute Russell's viper venom time for the diagnosis of lupus anticoagulant. Blood 1986; 68: 869-74.

Furie BC, Furie B.
Coagulant protein of Russell's viper venom. Methods Enzymol 1976; 45: 191-205.



Package size: Vial containing 5 U/vial
Vial containing 50 U/vial

Code: 121-06
121-07

FOR RESEARCH USE ONLY. NOT FOR HUMAN USE OR DRUG USE.

[Data sheet / 5802]
[RVV-X.doc]

®Registered Trademark of Pentapharm Ltd.

