



Product Data Sheet for:
Sheep anti-Hirudin
Purified IgG

Catalogue #: SAHD-IG
Lot #: XXXX
Expiry: XXXX

Immunogen: Purified recombinant full-length hirudin coupled to a carrier protein.

Presentation: XXXX ml of immunoglobulin purified from sheep antiserum. The IgG concentration is XXXX mg/ml in 10 mM HEPES, pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v). Total protein is **10 mg**.

Storage: Store between -10 and -20°C. Product will become viscous but will not freeze. Avoid storage in frost-free freezers. Keep vial tightly capped. Allow product to warm to room temperature and gently mix before use.

Specificity: This antibody is specific for hirudin as demonstrated by immunoelectrophoresis and ELISA.

Application: This product is suitable as a source of antibodies to hirudin.

For Research Use Only

Related products:

Product: SAHD-AP	Sheep anti-hirudin, affinity-purified IgG	0.5 mg vial
Product: SAHD-HRP	Sheep anti-hirudin, peroxidase conjugated purified IgG	0.2 mg vial

Visit our website (www.affinitybiologicals.com) for details.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Limited Warranty: This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.