



****REPRESENTATIVE DATASHEET****

Product Data Sheet for:
Sheep anti-human Neutrophil Elastase
Affinity-Purified IgG

Catalogue #: SANE-AP
Lot #: XXXX
Expiry: XXXX

Immunogen: Neutrophil elastase (active-site blocked) prepared from human sputum.

Presentation: XXXX ml of sheep immunoglobulin affinity purified on immobilized elastase (active site blocked) by ion exchange chromatography. The APIgG concentration is XXXX mg/ml in 10 mM HEPES, pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v). Total protein is **0.5 mg**.

Storage: Store between -10 and -20°C. Product will become viscous but will not freeze. Avoid storage in frost-free freezers. Keep vial tightly capped. Allow product to warm to room temperature and gently mix before use.

Specificity: This antibody is specific for neutrophil elastase as demonstrated by immunoelectrophoresis and ELISA.

Applications: This product is suitable as a source of enriched antibodies to neutrophil elastase.

For Research Use Only

Related products:

Product: SANE-IG	Sheep anti-human Neutrophil Elastase, purified IgG	10 mg vial
Product: SANE-HRP	Sheep anti-human Neutrophil Elastase, peroxidase conjugated IgG	0.2 mg vial

Visit our website (www.affinitybiologicals.com) for details.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Limited Warranty: This product is warranted to perform in accordance with its labelling and literature. Affinity Biologicals Inc. disclaims any implied warranty of merchantability or fitness for any other purposes, and in no event will Affinity Biologicals Inc. be liable for any consequential damages arising out of aforesaid express warranty.