



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Prekallikrein

Concentration: 1.96 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.510 ml

Activity: 28.00 PEU/mg

Total protein: 1.00 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M
NaCl/ pH 5.3

Storage: < -60 °C

Expiration Date: 02/2007

PO: N/A

The Human Prekallikrein is >95% pure as judged by SDS-PAGE gels and shows no reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma that was tested and found negative for HIV-1, HIV-2, Hepatitis B Surface antigen and HCV. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample